

BS

Surgical Instruments

DENTAL INSTRUMENTS CATALOGUE



www.bssurgicalinstruments.co.uk

about us.

At BS Surgical Instruments, our mission is to empower dental and surgical professionals with instruments that elevate patient care and clinical precision. We are proud to offer a comprehensive range of high-quality instruments designed to meet the day-to-day needs of modern hospitals and practices.

Our commitment to excellence is validated by our ISO and UK CA certifications, ensuring every instrument meets stringent international standards. Each product is developed under the careful watch of experienced experts and a strict quality control system, leaving no room for compromise.

Driven by a passion to make a positive impact, we continuously strive to innovate and enhance our service. At BS Surgical Instruments, every detail is focused on delivering value, supporting your clinical success, and ensuring that you can depend on our solutions every day. We express our heartfelt thanks to our customers for their unwavering trust and are dedicated to continuing to earn that confidence through exceptional service and unparalleled quality.

We believe that every tool we create contributes to a broader mission: improving healthcare outcomes and supporting the professionals who make it happen.


Surgical Instruments

Precision in Every Procedure




PROBES

HANDLE OPTIONS



EXPLORERS



HANDLES FOR MIRRORS

Handles with M2.5 Thread




MOUTH MIRRORS

Note: All Mouth Mirrors are Available in Plane & Magnifying Glass



Pack of 12



DRESSING PLIERS



CURETTES

GRACEY STANDARD

HANDLE OPTIONS



Silicone Rubber 9.5 mm



Colour Resin 9.5 mm



Stainless Steel 9.5 mm



BS-0134
Gracey 1-2



BS-0135
Gracey 7-8



BS-0136
Gracey 13-14



BS-0137
Gracey 3-4



BS-0138
Gracey 9-10



BS-0139
Gracey 15-16



BS-0140
Gracey 5-6



BS-0141
Gracey 11-12



BS-0142
Gracey 17-18



Gracey 5-6

**Application field of original
Gracey Curettes**

1-2, 3-4, 5-6
Front teeth

7-8, 9-10
Bicuspid and molars

11-12, 15-16
Bicuspid and molars mesial,
furcation: mesial
part of the distal root

13-14, 17-18
Bicuspid and molars distal,
furcation: distal part of the mesial
root.

UNIVERSAL CURETTES



8 mm



8 mm

HANDLE OPTIONS



BS-0143
Barnhart 1-2



BS-0144
McCall 13-14



BS-0145
McCall 17-18 S



BS-0146
Columbia 2L-2R



BS-0147
Langer 1-2



BS-0148
Barnhart 5-6



BS-0149
McCall 13-14 S



BS-0150
Columbia 4L-4R



BS-0151
Langer 3-4



BS-0152
Rules 3-4



BS-0153
McCall 17-18



BS-0154
Younger-Good 7-8



BS-0155
Langer 5-6



Columbia 4L-4R

SCALERS

HANDLE OPTIONS



6. mm



6 mm



9.5 mm



6.5 mm



BS-0156
Jacquette 3-4



BS-0159
Jacquette 14-15



BS-0162
Jacquette 3-U5



BS-0157
Jacquette 3-U15



BS-0160
Jacquette 4-U5



BS-0163
Jacquette 4-U15




BS-0158
M23



BS-0161
U15-30



BS-0164
U15-33



8 mm

H 6 -7



BS-0165
Eage Claw



BS-0170
Eagle Talon



BS-0166
Curaler - N 137



BS-0171
N135



BS-0167
S 311-312



BS-0172
N1



BS-0168
N3



BS-0173
SE Nebraska 128
(Single End)



BS-0169
Morse 0-00



BS-0174
N2

BS Surgical instrument is crafted by hand. Tips are fabricated individually by highly skilled precision technicians. Each instrument is continually held to strict quality controls during each stage of manufacturing, with tapering point working end. Well suitable for inter-dental cavities.

PLASTIC FILLING CARVERS



6 mm



9.5 mm



BS-0175
PF 11



BS-0179
PFI 8A



BS-0176
PFI ASH 6



BS-0180
PFI 13



BS-0177
PFI 6



BS-0181
PFI S6



BS-0178
PFI 11A



BS-0182
PFI D3

Original (rigid) with sharp cutting edges for carving. Our high grade steel and state-of-the art treatment processes provide highly polished surface that offer a non-stick surface.



5.5 mm PF 1 ASH 6

PLASTIC FILLING

TITANIUM NITRIDE COATING IS SIGNIFICANTLY HARDER THAN STEEL.



BS-0183
Composite 4
Paddle / Paddle




BS-0186
Composite Gregg 4-5
Paddle / Paddle




BS-0184
Composite 7
Paddle / Plugger

2.0 mm




BS-0187
Composite 6
Paddle / Plugger

1.0 mm



BS-0185
Composite 5

1.0 mm



BS-0188
Composite 6-21A
Paddle / Acorn

1.4 mm



0.9 mm

BS-0189 M-1 Ball Burnisher (Posterior)



BS-0190 Composite 9
Anterior






6 mm Stainless Steel Double End (Standard)




BS-0191
Wedelsteadt 3-4
(11.5-15-3)



BS-0192
Binangle 8-9
(20-9-8)




BS-0193
Wedelsteadt 5-6
(15-15-3)



BS-0194
Binangle 11-12
(15-8-8)

BS Surgical excavators, used for the removal of softened carious dentin, are equipped with spoon-shaped working end for work on deep areas and a round working end for tissue-saving cavity preparation.

EXCAVATORS



6 mm

1.2 mm



2.3 mm

BS-0195
17




1.5 mm

BS-0196
20



1.5 mm

BS-0197
42-43
(15-8-14)



1.2 mm

BS-0198
70-71



1.5 mm

BS-0199
E-1
(12-9-15)




1.15 mm

BS-0200
18



BS-0201
38-39
(11.5-7-14)



6.6 mm

1.0 mm



1.9 mm

BS-0202
64-65



BS-0203
74-75




1.5 mm

BS-0204
E-2
(15-9-15)



1.1 mm

BS-0205
40-41
(11-7-14)



1.0 mm

BS-0206
66-67




2.5 mm

BS-0207
78-79



2.5 mm

BS-0208
E-3
(25-9-15)




8 mm


74-75

CONDENSERS / PLUGGERS (NON-SERRATED)


HANDLE OPTIONS




BURNISHERS




CARVERS




KNIVES



WAXING INSTRUMENTS



MARGIN TRIMMERS



ENDODONTICS SPREADERS / PLUGGERS




CORD PACKERS

HANDLE OPTIONS






DISCOID-CLEOID




All of our instruments are hand sharpened by highly trained technicians to provide sharper, longer lasting tips and edges.

SINUS LIFT ELEVATORS



OSTEOTOMES

Handle Size: 10 mm, 13 mm



CHEEK, LIP RETRACTOR AND STANDS



PERIOSTEALS ELEVATOR


For reflecting and retracting the mucoperiosteum after incisions of the gingival tissue.



BONE SPLITTING INSTRUMENTS



BONE FILES





BS-321
Miller & Lucas
Spoon Shape
2.2mm, 2.8 mm
6

BS-322
Will
2.5mm



BS-313
Articulating Paper Forcep
Miller
15 cm



BS-314
Williger 14.5 cm

Sharp Spoon



BS-315
Hemingwa
16.5 cm



BS-316
Will Sharp Spoon 16.5 cm

2,5 mm

3,5 mm



BS-317
Periosteal Elevator Freer (Sharp / Blunt)



BS-318
Bone Chisel Partsch 17 cm (3-6 mm)



Partsch **BS-319**
Gouge Width 4 mm, 5 mm, 6 mm, 7 mm

17 cm (3-6 mm)



BS-320
Gouge
17 cm (3-6 mm)

WATER AND CARTRIDGE SYRINGES

Also Available With Aspiration



BS-0332
Aspirating Cartridge Syringe
metric / with US thread
screw tops for cannulas ointment cannula
blunt/curved
1.8 cc, 2.2 cc



BS-0333
Water Syringe
10 ml with Luer
lock cannula



BS-0334
Cartridge Syringe
(TKC), cylinder
without aspiration metric / US
threads chrome plated brass
1.8 cc, 2.2 cc



BS-0335
Skypject Intra lingual Syringe
Used in performing PDL
injections. Brass
and stainless steel
11 cm



BS-0336
Self Aspirating Syringe
anodized aluminum,
and stainless steel
1.8 cc, 2.2 cc

MEASUREMENT / ENDODONTIC

MEASURING INSTRUMENTS



BS-0337
Positive Tube Holder
Anodized
Aluminium



BS-0338
Boley Caliper
with Lock 14 cm
Scale 0-100 mm



BS-0339
Measurement of
bone thickness with a precision
of up to 0.05 mm
14.5 cm



BS-0340
Ivanson
Wax/Metal 10 cm (Scale 0-10)
New Dental Caliper



BS-0341
Castroviejo Caliper
10 cm (Scale 0-40)

RUBBER DAM INSTRUMENTS



BS-0342
Rubber dam punch
Anisworth 16 cm
also available Ivory with 1 and 2 rows of holes




BS-0343
Rubber dam
clamp forcep Brewer
17 cm



BS-0344
Rubber dam clamp forcep
Stokes
17 cm

RUBBER DAM CLAMPS **BS-C345**



BS-F345-A
Rubber dam frame
big/small 15X15 cm
Stainless Steel/Plastic

Assorted Colours Moulded Plastic


Pack of 6



BS-0356



BS-0357



BS-0358

		mm	mm
U0	xxs	63	48
U1	xs	65	52
U2	s	68	54
U3	M	71	58
U4	L	73	59
U5	XL	78	63

Perforated Without Retention Rim



BS-0359

BS-0360

Perforated With Retention Rim



BS-0361

BS-0362

		mm	mm
L0	xxs	56	44
L1	xs	66	50
L2	s	70	52
L3	M	73	56
L4	L	76	58
L5	XL	78	60

Un perforated With Retention Rim



BS-0363

BS-0364

Also Available Partial Trays



BS-0365
Perforated for
edentulous upper jaws

BS-0366
Un perforated For
toothed lower jaws

		mm	mm
L		79	70
M		76	67
S		73	60

		mm	mm
L		79	70
M		76	67
S		73	60



BS-0367
Perforated
For Partial
Impressions right
R1, R2, R3



BS-0368
Unperforated
For Partial
Impressions Left
L1, L2, L3

AMALGUM CARRIER, WELL AND GUN



BS-0369
23 cm

BS-0370
Ivory 23 cm
Also Available Single Ended 17.5 cm



BS-0371
Amalgum Well
H 1.5 cm / D 3.5 cm
H 3.5 cm / D 3.5 cm



13 cm, 15 cm
BS-0372
Amalgum Gun
30° 45° 90°

ENGLISH PATTERN



BS-0373
Fig. 1
For Upper Incisors
and Canines, Wide



BS-0374
Fig. 168
For Upper Canines



BS-0375
Fig. 136
For Upper Premolars



BS-0376
Fig. 17
For Upper Molars, Right



BS-0377
Fig. 18
For Upper Molars, Left



BS-0378
Fig. 7
For Upper Premolars



BS-0379
Fig. 67
For Upper Wisdom Teeth, Narrow
Available: For Upper Wisdom Teeth, Wide



BS-0380
Fig. 29
For Upper Front Roots



BS-0381
Fig. 30
For Upper Roots



FORCEPS FOR CHILDREN



BS-0385
Fig. 89
For Upper Molars,
Right, Gripping In Depth



BS-0382
Fig. 18 A
For Upper Molars, Either Side



BS-0383
Fig. 51 A
For Upper Front Roots,
Narrow Beaks



BS-0384
Fig. 13
For Lower Canines And Premolars



BS-0386
Fig. 33 A
For Lower Roots, With
Narrow Beaks



BS-0387
Fig. 33
For Lower Incisors,
Premolars And Roots



BS-0388
Fig. 75
For Lower Canines
And Premolars



BS-0389
Fig. 24
For Lower Molars, Left



BS-0390
Fig. 31
For Lower Premolars
And Roots



BS-0391
Fig. 73
For Lower Molars



BS-0392
Fig. 21
For Lower Wisdom
Teeth And Molars

BSS

Surgical Instruments



**Scan for
PDF Catalogue**



**Scan for
website**

Surgical Instruments



Unit 1 Vulcan centre Dixon st WV2 2BS
Wolverhampton



004420 3889 4690



sales@bssurgicalinstruments.co.uk



www.bssurgicalinstruments.co.uk